

www.ixcellsbiotech.com

10340 Camino Santa Fe, Suite C, San Diego, CA 92121 Tel: (858)412-5988 Fax: (858)368-8716 Technical Supports: <u>supports@ixcellsbiotech.com</u> Orders: <u>orders@ixcellsbiotech.com</u>

# **Product Information**

### Human Brain Microvascular Endothelial Cells

Catalog Number	10HU-051	Cell Number	0.5 x 10 <sup>6</sup> cells/vial
Species	Homo sapiens	Storage Temperature	Liquid Nitrogen

# Description

The human brain microvascular endothelial cells (HBMVEC) are the major element of the blood-brain barrier and comprise the primary limitation to passage of substances, both soluble and cellular, from the blood into the brain. HBMVEC utilizes unique features that distinguish themselves from those of peripheral endothelial cells. Most prominent among these are the following: 1) intercellular 'tight junctions' that display high transendothelial electrical resistance and retard paracellular flux [1], 2) absence of fenestrae and a reduced level of fluid-phase endocytosis [2], and 3) asymmetrically-localized enzymes and carrier-mediated transport systems that engender a truly 'polarized' phenotype. Like peripheral endothelial cells, however, HBMVEC express, or can be induced to express, cell adhesion molecules on their surface that regulate the extravasation of leukocytes into the brain. HBMVEC has been widely used for studying the molecular and cellular basis of blood-brain barrier.

iXCells Biotechnologies provides high quality HBMVEC, which are isolated from human brain and cryopreserved at P1, with >0.5 million cells in each vial. These HBMVEC express vWF/Factor VIII, CD31 (PECAM), and Dil-Ac-LDL by uptake. They are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast, and fungi and can further expand for 10 population doublings in Endothelial Cell Growth Medium (Cat# MD-0010) under the condition suggested by iXCells Biotechnologies.



Figure 1. Human Brain Microvascular Endothelial Cells (HBMVEC), phase contrast.

For Research Only

1

All Rights Reserved

iXCells Biotechnologies USA, LLC.



Figure 2. Immunostaining of Human Brain Microvascular Endothelial Cells (HBMVEC) using the antibodies against vWF (green) and CD31 (red).

# **Product Details**

Tissue	Human brain
Package Size	0.5 x 10 <sup>6</sup> cells/vial
Passage Number	P1
Shipped	Cryopreserved
Storage	Liquid nitrogen
Growth Properties	Adherent
Media	Endothelial Cell Growth Medium (Cat# MD-0010)

## **Protocols**

2

### **Thawing of Frozen Cells**

- 1. Upon receipt of the frozen cells, it is recommended to thaw the cells and initiate the culture immediately in order to retain the highest cell viability.
- 2. To thaw the cells, put the vial in 37°C water bath with gentle agitation for ~1 minute. Keep the cap out of water to minimize the risk of contamination.
- 3. Pipette the cells into a 15ml conical tube with 5ml fresh Endothelial Cell Growth Medium (Cat# MD-0010).

All Rights Reserved

- 4. Centrifuge at 1000rpm (~220g) for 5 minutes under room temperature.
- 5. Remove the supernatant and resuspend the cells in Endothelial Cell Growth Medium (Cat# MD-0010).
- 6. Culture the cell in T75 flask.

Safety Precaution: it is highly recommended that protective gloves and clothing should be used when handling frozen vials.

#### **Standard Culture Procedure**

- **1.** HBMVEC can be cultured in Endothelial Cell Growth Medium (Cat# MD-0010).
- 2. When cells reach ~80-90% confluence, remove the medium, and wash once with sterile PBS (5ml/T75 flask).
- Add ~2.5ml of 0.05% Trypsin-EDTA to the flask and incubate for ~3 minutes at 37°C. Neutralize the enzyme by adding 2-3 volumes of cell culture medium.
- 4. Centrifuge 1,000rpm (~220g) for 5min and resuspend the cells in desired volume of medium.
- 5. Seed new culture vessels at  $5 \times 10^3$  cells/cm<sup>2</sup>.

### Reference

[1] Crone, C. and Oleson, S. P. (1992) Electrical resistance of brain microvessel endothelium. Brain Res. 241: 49-55.

[2] Reese, T. S. and Karnovsky, M. J. (1967) Fine structural localization of blood-brain barrier to exogenous peroxidase. J. Cell Biol. 34:9-14. [3] Vorbrodt, A. W. (1988) Ultrastructural cytochemistry of blood-brain barrier endothelia. Prog. Histochem. Cytochem. 18(3):1-96.

#### **Disclaimers**

3

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures. © iXCells Biotechnologies 2015. All rights reserved.

For Research Only

All Rights Reserved

iXCells Biotechnologies USA, LLC.