

Product Information

Human Epidermal Keratinocyte-adult (HEK-a)

Catalog Number	10HU-225	Cell Number	0.5 x 10 ⁶ cells/vial
Species	<i>Homo sapiens</i>	Storage Temperature	Liquid Nitrogen

Description

A keratinocyte is the predominant cell type in the epidermis, constituting 90% of the cells found there^[1]. The primary function of keratinocytes is the formation of a barrier against environmental damage by pathogenic bacteria, fungi, parasites, viruses, heat, UV radiation and water loss. Keratinocytes are also able to produce a variety of cytokines, growth factors, interleukins and complement factors. Therefore keratinocytes are important for wound healing, inflammation, and immune response.

iXCells Biotechnologies provides high quality Human Epidermal Keratinocyte-adult (HEK-a), which are isolated from adult skin samples and cryopreserved at P1, with >0.5 million cells in each vial. HEK-a are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast, and fungi. They can further expand for 10 population doublings in Keratinocyte Growth Medium (Cat# MD-0047) under the condition suggested by iXCells Biotechnologies.

Product Details

Tissue	Adult skin
Package Size	0.5x10 ⁶ cells/vial
Passage Number	P1
Shipped	Cryopreserved
Storage	Liquid nitrogen
Growth Properties	Adherent
Media	Keratinocyte Growth Medium (MD-0047)

Protocols

Thawing of Frozen Cells

1. Upon receipt of the frozen cells, it is recommended to thaw the cells and initiate the culture immediately in order to retain the highest cell viability.
2. To thaw the cells, put the vial in 37°C water bath with gentle agitation for ~1 minute. Keep the cap out of water to minimize the risk of contamination.
3. Pipette the cells into a 15 ml conical tube with 5ml fresh Keratinocyte Growth Medium (MD-0047).
4. Centrifuge at 1,000 rpm (~220 g) for 5 minutes under room temperature.
5. Remove the supernatant and resuspend the cells in fresh culture medium.
6. Culture the cell in T75 flask.

Safety Precaution: *it is highly recommended that protective gloves and clothing should be used when handling frozen vials.*

Standard Culture Procedure

1. Keratinocytes can be cultured in Keratinocyte Growth Medium (MD-0047).
2. When cells reach ~80-90% confluence, remove the medium, and wash once with sterile PBS (5ml/T75 flask).
3. Add ~2.5ml of 0.25% Trypsin-EDTA to the flask and incubate for ~3 minutes at 37°C. Neutralize by adding 2-3 volumes of Keratinocyte Growth Medium.
4. Centrifuge 1,000 rpm (~220 g) for 5min and resuspend the cells in desired volume of Keratinocyte Growth Medium.
5. Seed the cells onto new culture vessels at 3×10^3 cells/cm².

Reference

[1] McGrath JA, Eady RAJ, Pope FM. (2004). "Anatomy and Organization of Human Skin". In Burns T, Breathnach S, Cox N, Griffiths C. Rook's Textbook of Dermatology (7th ed.). Blackwell Publishing. p. 4190.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures.

© iXCells Biotechnologies 2015. All rights reserved.