

10340 Camino Santa Fe, Suite C, San Diego, CA 92121 Tel: (858)412-5988 Fax: (858)368-8716

Technical Supports: <a href="mailto:supports@ixcellsbiotech.com">supports@ixcellsbiotech.com</a>
Orders: orders@ixcellsbiotech.com

# **Product Information**

### **Human Intestinal Microvascular Endothelial Cells (HIMEC)**

Catalog Number	10HU-065	Cell Number	0.5 million cells/vial
Species	Homo sapiens	Storage Temperature	Liquid Nitrogen

# **Description**

Endothelial cells lining the microvasculature are known to play a critical "gatekeeper" role in the inflammatory process through their ability to recruit circulating immune cells into tissues and foci of inflammation. Studies show that intestinal microvascular endothelial cells (IMEC) exhibit a strong immune response to LPS challenge and play a critical regulatory role in gut inflammation [1]. Pharmacological inhibition of NOS in activated HIMEC resulted in a significant increase in leukocyte binding [2]. Gene expression profile studies revealed that intestinal endothelial cells express biotinidase, which is involved in biotin recycling [3]. HIMEC cultures have enabled scientists to perform systematic analyses of cytokine profiles with regard to mRNA expression and protein secretion, and to compare these data with cytokine profiles from other endothelial cells.

iXCells Biotechnologies provides high quality HIMEC, which are isolated from human intestinal tissue and cryopreserved at P1, with >0.5 million cells in each vial. HIMEC express vWF/Factor VIII, CD31 (PECAM), and Dil-Ac-LDL by uptake. They are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast, and fungi and can further expand in Endothelial Cell Growth Medium (Cat# MD-0010) under the condition suggested by iXCells Biotechnologies.

## **Product Details**

Tissue	Normal human intestine tissue	
Package Size	0.5 million cells/vial	
Passage Number	P1	
Shipped	Cryopreserved	
Storage	Liquid nitrogen	
<b>Growth Properties</b>	Adherent	
Media	Endothelial Cell Growth Medium (Cat# MD-0010)	

## **Protocols**

### **Thawing of Frozen Cells**

- Upon receipt of the frozen cells, it is recommended to thaw the cells and initiate the culture immediately in order to retain the highest cell viability.
- 2. To thaw the cells, put the vial in 37°C water bath with gentle agitation for 1-2 minutes. Keep the cap out of water to minimize the risk of contamination.
- 3. Pipette the cells into a 15 mL conical tube with 5ml fresh Endothelial Cell Growth Medium (Cat# MD-0010).
- 4. Centrifuge at 1,000 rpm (~220 g) for 5 minutes under room temperature.
- 5. Remove the supernatant and resuspend the cells in fresh Endothelial Cell Growth Medium.
- 6. Culture the cell in the T75 flask. Change the medium every other day until cells reach 80-90% confluence.

Safety Precaution: it is highly recommended that protective gloves and clothing should be used when handling frozen vials.

#### Standard Culture Procedure

- 1. HIMEC can be cultured in Endothelial Cell Growth Medium (Cat# MD-0010).
- 2. When cells reach ~80-90% confluence, remove the medium, and wash once with sterile PBS (5mL for one T75 flask).
- 3. Add 3 mL of 0.25% Trypsin-EDTA to the flask and incubate for 5 minutes at 37°C. Neutralize the enzyme by adding 2-3 volumes of cell culture medium.
- 4. Centrifuge 1,000 rpm (~220 g) for 5 minutes and resuspend the cells in desired volume of medium.
- Seed the cells in the new culture vessels at 5 x 10<sup>3</sup> cells/cm<sup>2</sup>. Change the medium every other day until cells reach 80-90% confluence.

#### **Disclaimers**

This product is intended for laboratory research purposes only. It is not intended for use in humans. While iXCells Biotechnologies uses reasonable efforts to include accurate and up-to-date information on this product sheet, we make no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. iXCells Biotechnologies does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, and use. iXCells Biotechnologies is not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, iXCells Biotechnologies is not liable for damages arising from the misidentification or misrepresentation of cultures.

© iXCells Biotechnologies 2015. All rights reserved.